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nature biotechnology

Explore content ▾ About the journal ▾ Publish with us ▾ Subscribe

[nature](#) > [nature biotechnology](#) > [letters](#) > article

Letter | [Published: 11 May 2020](#)

Multiplex digital spatial profiling of proteins and RNA in fixed tissue

[Christopher R. Merritt](#) , [Giang T. Ong](#), [Sarah E. Church](#), [Kristi Barker](#), [Patrick Danaher](#), [Gary Geiss](#), [Margaret Hoang](#), [Jaemyeong Jung](#), [Yan Liang](#), [Jill McKay-Fleisch](#), [Karen Nguyen](#), [Zach Norgaard](#), [Kristina Sorg](#), [Isaac Sprague](#), [Charles Warren](#), [Sarah Warren](#), [Philippa J. Webster](#), [Zoey Zhou](#), [Daniel R. Zollinger](#), [Dwayne L. Dunaway](#), [Gordon B. Mills](#) & [Joseph M. Beechem](#) 

Nature Biotechnology **38**, 586–599 (2020) | [Cite this article](#)

20k Accesses | 119 Citations | 133 Altmetric | [Metrics](#)

Abstract

Digital Spatial Profiling (DSP) is a method for highly multiplex spatial profiling of proteins or RNAs suitable for use on formalin-fixed, paraffin-embedded (FFPE) samples. The approach relies on (1) multiplexed readout of proteins or RNAs using oligonucleotide tags; (2) oligonucleotide tags attached to affinity reagents (antibodies or RNA probes) through a photocleavable (PC) linker; and (3) photocleaving light projected onto the tissue sample to release PC oligonucleotides in any spatial pattern across a region of interest (ROI) covering 1 to ~5,000 cells. DSP is capable of single-cell sensitivity within an ROI using the antibody readout, with RNA detection feasible down to ~600 individual mRNA transcripts. We show spatial profiling of up to 44 proteins and 96 genes (928 RNA probes) in lymphoid, colorectal tumor and autoimmune tissues by using the nCounter system and 1,412 genes (4,998 RNA probes) by using next-generation sequencing (NGS). DSP may be used to profile not only proteins and RNAs in biobanked samples but also immune markers in patient samples, with potential prognostic and predictive potential for clinical decision-making.

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Implementation of a Multiplex and Quantitative Proteomics Platform for Assessing Protein Lysates Using DNA-Barcoded Antibodies[†]

Jinho Lee¹, Gary K. Geiss¹, Gokhan Demirkan¹, Christopher P. Vellano¹, Brian Filanoski¹, Yiling Lu¹, Zheslin Ju¹, Shuangxing Yui¹, Huihang Guo¹, Lisa Y. Bogatzki¹, Warren Carter¹, Rhonda K. Meredith¹, Savitri Krishnamurthy¹, Zhiyong Dong¹, Joseph M. Beechem¹, and Gordon B. Mills¹

Multiplex digital spatial profiling of proteins and RNA in fixed tissue

Christopher B. Moritt¹, Giang T. Ong¹, Sarah E. Church¹, Kristi Barker¹, Patrick Danaher¹, Gary Geiss¹, Margaret Hoang¹, Jaemyeong Jung¹, Yan Liang¹, Jill McKay-Fleisch¹, Karen Nguyen¹, Zach Norgaard¹, Kristina Sorg¹, Isaac Sprague¹, Charles Warren¹, Sarah Warren¹, Philippe J. Webster¹, Zoey Zhou¹, Daniel R. Zollinger¹, Dwayne L. Dunaway¹, Gordon B. Mills¹ and Joseph M. Beechem^{1,2}

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

GEP:
Multi-plex PCR (bulk),
RNAseq (bulk),
Single-cell RNAseq,
PEP:
Flow Cytometry,
Mass cytometry,
DSP

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CLINICAL CANCER RESEARCH

High-Plex Predictive Marker Discovery for Melanoma Immunotherapy-Treated Patients Using Digital Spatial Profiling

CLINICAL CANCER RESEARCH

Biomarkers Associated with Beneficial PD-1 Checkpoint Blockade in Non-Small Cell Lung Cancer (NSCLC) Identified Using High-Plex Digital Spatial Profiling

nature cancer

Spatial proteomic characterization of HER2-positive breast tumors through neoadjuvant therapy predicts response

ESMO

High multiplex analysis of the immune microenvironment in bone marrow trephine samples using GeoMx™ digital spatial profiling

cancers


High-Plex and High-Throughput Digital Spatial Profiling of Non-Small-Cell Lung Cancer (NSCLC)

James Monkman, Tomraj Tabiri, Majid Ebrahimi Warkiani, Connor O'Leary, Rabal Ludwa, Derek Richard, Ken O'Byrne, and Aratha Kalasinghe

IOTECH

High multiplex analysis of the immune microenvironment in bone marrow trephine samples using GeoMx™ digital spatial profiling

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Journal for Immunotherapy of Cancer

Distinct tumor microenvironments of lytic and blastic bone metastases in prostate cancer patients

Qian L. He, Meredith D. Prosser, Christine M. Seignin, E. Din Smith, Susan M. Edgermon, Adrie Van Balkhoven, M. Scott Lucia, and Philip Ouyang

Immunologic alterations in the pancreatic cancer microenvironment of patients treated with neoadjuvant chemotherapy and radiotherapy

Matthew R. Farnen, Layal Sayegh, Michael Brandon Ware, Hsiao-Rong Chen, Jingling Gong, Yan Liang, Alyssa Kravitskaya, Shichik K. Mutha, Mohammad Zaidi, Juan M. Samir, David Gaddy, Patrick Patel, Bassem El-Rayes, Walid Ibrahim, and Gregory B. Leshem

Journal for Immunotherapy of Cancer

Inter- and intra-tumor heterogeneity of metastatic prostate cancer determined by digital spatial gene expression profiling

Lauren Brady, Michelle Kriner, Isaac Coleman, Cole Morrissey, Martine Roudier, Lawrence D. True, Roman Gulati, Stephen R. Plymate, Zoey Zhou, Brian Biditt, Rhonda Meredith, Gary Geiss, Margaret Hoang, Joseph Beechem, and Peter S. Nelson

High dose-rate brachytherapy of localized prostate cancer converts tumors from cold to hot

Simon P. Kaam, Heiwoo Hahn, Thi Nguyen, Minyu Wang, Nicholas Van Kesteren, Catherine Mitchell, Franco Carraia, David J Byrne, Sue Haupt, Georgina Flynn, Phillip K Darcy, Shaheen Sandhu, Piers Blombery, Ygal Haupt, Scott G Williams, Paul J Neeson

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nature medicine

Opposing immune and genetic mechanisms shape oncogenic programs in synovial sarcoma

Livnat Jerby-Aron, Cyril Nefzi, Marri E. Shon, Hannah R. Weisman, Nathan D. Mathewson, Matthew J. McBride, Brian Hays, Benjamin Izan, Angela Volker, Gayle Boulay, Luisa Ghent, Alyssa R. Richman, Liliane C. Breyer

Onco-fetal Reprogramming of Endothelial Cells Drives Immunosuppressive Macrophages in Hepatocellular Carcinoma

Ankur Sharma, Justine Jia Wen Seow

Cell

Spatially organized multicellular immune hubs in human colorectal cancer

scRNA-seq and spatial profiling of CRC patients

Tertiary lymphoid structures improve immunotherapy and survival in melanoma

Wito Cabrita, Martin Lasser, Adriana Serrao, Marco Donati, Mathilde Stangor, Shreshth Mittal, Ivo Johansson, Bengt Philip, Kristina Herberich, John Wilton Christensen, Alben van Schrick, Katharina Langner, Sarah Warner, Karim Alnaitan, Mikael Olsson, Kristian Preker, Christian Ingemar, Karoline Sultmann, Erik Schindler, Henrik Schmalz, Lars Bastholt, Ana Carmona, Jennifer A. Wargo, Inga Marie Svane, and Ola Jonsson

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CLINICAL CANCER RESEARCH
High-Plex Predictive Marker Discovery for Melanoma Immunotherapy-Treated Patients Using Digital Spatial Profiling

nature cancer
Spatial proteomic characterization of HER2-positive breast tumors through neoadjuvant therapy predicts response

cancers
High-Plex and High-Throughput Digital Spatial Profiling of Non-Small-Cell Lung Cancer (NSCLC)

CLINICAL CANCER RESEARCH
Biomarkers Associated with Beneficial PD-1 Checkpoint Blockade in Non-Small Cell Lung Cancer (NSCLC) Identified Using High-Plex Digital Spatial Profiling


ESMO
High multiplex analysis of the immune microenvironment in bone marrow trephine samples using GeoMx™ digital spatial profiling


IOTECH
High multiplex analysis of the immune microenvironment in bone marrow trephine samples using GeoMx™ digital spatial profiling

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

Journal for Immunotherapy of Cancer
Distinct tumor microenvironments of lytic and blastic bone metastases in prostate cancer patients

ARTICLE
Inter- and intra-tumor heterogeneity of metastatic prostate cancer determined by digital spatial gene expression profiling

JCI INSIGHT
Immunologic alterations in the pancreatic cancer microenvironment of patients treated with neoadjuvant chemotherapy and radiotherapy


Open access
High dose-rate brachytherapy of localized prostate cancer converts tumors from cold to hot

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

nature medicine
Opposing immune and genetic mechanisms shape oncogenic programs in synovial sarcoma

Cell
Spatially organized multicellular immune hubs in human colorectal cancer

Cell
Onco-fetal Reprogramming of Endothelial Cells Drives Immunosuppressive Macrophages in Hepatocellular Carcinoma

Article
Tertiary lymphoid structures improve immunotherapy and survival in melanoma

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




Journal of Immunotherapy of Cancer 2022 (肉瘤)

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Table 1 Summary of patient and tumor characteristics

| Patient | Tissue | Diagnosis | ER(%)/ PR(%)/ HER2(+) | PAM50 Subtype | TIS Status | Treatment at the time of sample collection |
|---------|-------------|------------------|-----------------------------|---------------|------------|--|
| 1 | Breast | IDC, DCIS | 2/0/3 | Luminal A | High | Gemcitabine and trastuzumab |
| 1 | Breast | Tumor emboli | N/A | Luminal A | High | Nab-paclitaxel and trastuzumab |
| 1 | Lung | Metastatic | 0/0/3 | HER2 enriched | High | Nab-paclitaxel and trastuzumab |
| 2 | Breast | ICD, DCIS | 0/0/3 | HER2 enriched | Low | None. Initial diagnosis |
| 2 | Breast | Residual IDC | 0/0/N/A | HER2 enriched | High | TDM1, pertuzumab, doxorubicin and cyclophosphamide* |
| 2 | Soft tissue | Local recurrence | 0/0/3 | HER2 enriched | Low | Trastuzumab and pertuzumab |
| 3 | Brain | Metastatic | 15/0/3 | HER2 enriched | Low | None. Discontinued due to toxicities |
| 4 | Breast | IDC | 0/0/3 | Basal like | High | None. Initial diagnosis |
| 4 | Breast | Local recurrence | N/A | Luminal A | High | Doxorubicin, cyclophosphamide, paclitaxel and trastuzumab |
| 5 | Brain | Metastatic | 0/0/3 | Basal like | Low | TDM1 and anastrozole |
| 6 | Breast | IDC, DCIS | 80/70/2** | Luminal A | High | None. Initial diagnosis |
| 7 | Breast | IDC, DCIS | 0/0/3 | HER2 enriched | High | None. Initial diagnosis |
| 7 | Breast | Local recurrence | N/A | Basal like | High | None. Completed neoadjuvant and adjuvant therapy |
| 8 | Breast | IDC, DCIS | 10/15/1 | Luminal A | Low | None. Initial diagnosis |
| 8 | Soft tissue | Local recurrence | 0/0/3 | Basal like | Low | Nab-paclitaxel and atezolizumab (previously had triple negative disease) |

DCIS: ductal carcinoma in situ; IDC: ductal carcinoma in situ; ER: estrogen receptor; HER2: human epidermal growth factor 2; IDC: invasive ductal carcinoma; IDC: invasive ductal carcinoma; PR: progesterone receptor; TDM1: ado-trastuzumab emtansine; TIS: tumor inflammation score

* Neoadjuvant therapy in the ISPY2 clinical trial

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






Five Unique Fluorescent Profiling Modalities

Geometric

CD3 PanCK DNA


Assess heterogeneity
(e.g. tumor margin vs. core, cortical layers)

Segmentation


PanCK DNA


Evaluate biological compartments
(e.g. Tumor-TME)?

Cell Type/Rare cell


Profile cell populations
(e.g. CD3+ cells, GFAP+ cells)

Contour

Determine proximal response
(e.g. CART, AD plaque)

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








Spatially-distinct changes in immune cell activation

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| Timepoint | Segment | PAM50 | ROI.Type | Tissue | TIS.Status |
|-----------|-----------|------------|----------|--------|------------|
| Primary 1 | PanCK neg | Basal like | Cold | Brain | High |
| Primary 2 | PanCK pos | HER2e | Warm | Breast | Low |
| Met | | Luminal A | Hot | Lung | |

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Visualized Cancer Medicine 2021, 2, 1
© The Authors, published by EDP Sciences, 2021
<https://doi.org/10.1051/vcm/2020002>

Available online at:
<https://www.edpsciences.org>

TECHNOLOGICAL ADVANCES OPEN ACCESS

Spatially-resolved proteomics and transcriptomics: An emerging digital spatial profiling approach for tumor microenvironment

Nan Wang¹, Rongshui Wang¹, Xue Zhang¹, Xia Li¹, Yan Liang², and Zhiyong Ding^{1,*}

¹ Mills Institute for Personalized Cancer Care, Fynn Biotechnologies Ltd., Gangxing 3rd Rd, High-Tech and Innovation Zone, Bldg. 2, Rm. 2201, Jinan City, Shandong Province 250101, PR China
² NanoString Technologies, Inc., Seattle, WA, USA

Received 17 September 2020, Accepted 4 December 2020, Published online 3 March 2021

Received: 5 February 2021 | Revised: 28 May 2021 | Accepted: 3 June 2021
DOI: 10.1002/iot.202100041

Biotechnology Journal

MINI-REVIEW

Spatial transcriptomics and proteomics technologies for deconvoluting the tumor microenvironment

Ding et al. *Molecular Cancer* (2022) 21:53
<https://doi.org/10.1186/s12943-022-01526-8>

Molecular Cancer

REVIEW Open Access

Proteomics technologies for cancer liquid biopsies

Zhiyong Ding^{1*}, Nan Wang¹, Ning Ji^{2,3} and Zhe-Sheng Chen^{2*}

Abstract

Journal of Ethnopharmacology 289 (2022) 115018

Contents lists available at ScienceDirect

Journal of Ethnopharmacology

journal homepage: www.elsevier.com/locate/jethpharm

Ameliorative effect of *Gastrodia elata* Blume extracts on depression in zebrafish and cellular models through modulating reticulum 4 receptors and apoptosis

Rongchun Wang^{b,c,e,1}, Qingyu Ren^{a,b,c,e,1}, Daili Gao^{b,c,e,1}, Yam Nath Paudel^d, Xia Li^e, Lizhen Wang^{b,c}, Pengyu Zhang^{a,b,c}, Baokun Wang^{b,c}, Xueliang Shang^{a,c}, Meng Jin^{b,c,e,*}

^a School of Psychology, North China University of Science and Technology, 21 Bohai Road, Tangshan, 063210, Hebei Province, PR China
^b Biology Institute, Qilu University of Technology (Shandong Academy of Sciences), 26789 East Jingchi Road, Jinan, 250103, Shandong Province, PR China
^c Key Laboratory for Drug Screening Technology (Shandong Academy of Sciences), 26789 East Jingchi Road, Jinan, 250103, Shandong Province, PR China
^d Neuropharmacology Research Strength, Jeffrey Cheah School of Medicine and Health Sciences, Monash University Malaysia, Bandar Sunway, 47500, Selangor, Malaysia
^e Mills Institute for Personalized Cancer Care, Fynn Biotechnologies Ltd., Gangxing 3rd Rd, High-Tech and Innovation Zone, Bldg. 2, Rm. 2201, Jinan, 250101, Shandong Province, PR China

Contents lists available at ScienceDirect

Molecular and Cellular Endocrinology

journal homepage: www.elsevier.com/locate/mce

EVI1 overexpression promotes ovarian cancer progression by regulating estrogen signaling

Zixiang Wang^{a,b}, Yingwei Li^a, Nan Wang^c, Peng Li^{b,**}, Beihua Kong^{b,***}, Zhaojian Liu^{b,*}

^a Department of Obstetrics and Gynecology, Qilu Hospital, Cheeoo College of Medicine, Shandong University, Jinan, Shandong, 250012, China
^b Department of Cell Biology, School of Basic Medical Sciences, Cheeoo College of Medicine, Shandong University, Jinan, 250012, China
^c Mills Institute for Personalized Cancer Care and Fynn Biotechnologies Ltd., Jinan, Shandong, 250012, China

Drug Discovery Today • Volume xxx, Number xx • xxxx 2021

REVIEW

Recent progress on targeting leukemia stem cells

Xiang-Yu Ma^a, Liuya Wei^{a,c}, Zining Lei^b, Yanglu Chen^c, Zhiyong Ding^d, Zhe-Sheng Chen^{a,c}

^a School of Pharmacy, Welfang Medical University, Welfang 261053, PR China
^b Department of Pharmaceutical Sciences, College of Pharmacy and Health Sciences, St. John's University, Queens, NY 11439, USA
^c Columbia University Vagelos College of Physicians and Surgeons, New York, NY 10032, USA
^d Mills Institute for Personalized Cancer Care, Fynn Biotechnologies Ltd., Gangxing 3rd Rd, High-Tech and Innovation Zone, Jinan, Shandong 250101, PR China

frontiers in Bioengineering and Biotechnology

BRIEF RESEARCH REPORT
published: 28 October 2021
doi: 10.3389/fbioe.2021.757378

Tumor Microenvironment Profiles Reveal Distinct Therapy-Oriented Proteogenomic Characteristics in Colorectal Cancer

Nan Wang¹, Rongshui Wang¹, Xia Li¹, Zhenhao Song¹, Lingbo Xia¹, Jue Wang², Li Zhang³, Aiwu Wu^{4*} and Zhiyong Ding^{1*}


¹ Mills Institute for Personalized Cancer Care, Fynn Biotechnologies Ltd., Jinan, China, ² Department of Oncology, Shandong Provincial Hospital Affiliated to Shandong First Medical University, Jinan, China, ³ Department of Pathology, Beijing Cancer Hospital, Beijing, China, ⁴ Key Laboratory of Carcinogenesis and Translational Research (Ministry of Education), Department of Gastrointestinal Surgery, Peking University, Beijing, China

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